



	Mask to save
1	DAPI/nucleus mask
2	Voronoi mask
3	Stroma mask
3	Voronoi with stroma subtraction
4	Cytoplasm mas (Voronoi after nucleus and stroma subtraction)
5	DAPI/nucleus mask w/ stroma subtraction



DAPI




Voronoi



Cell membrane
 Voronoi with stroma subtraction and nucleus subtraction



Stroma mask



Voronoi/Cell mask

Why do we quantify MenalNV in stoma only?

We exclude stroma because we want to assess MenalNV intensity in tumor cells only and the fibroblasts in stroma can express quite a lot of MenalNV.